


Product Datasheet

Angiotensinogen antibody orb357928

| | |
|----------------------------|--|
| Description | Rabbit polyclonal antibody to Angiotensinogen |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse |
| Conjugation | Unconjugated |
| Custom conjugations | Unconjugated, Biotin, FITC, HRP, CF647, CF488A, CF594, CF405M |
| Tested Applications | ELISA, IHC-P, WB |
| Immunogen | Recombinant human Angiotensinogen protein (35-285aa) |
| Target | Angiotensinogen |
| Form/Appearance | Liquid: PBS with 0.02% sodium azide, 50% glycerol, pH7.3. |
| Storage | Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles. |
| Note | For research use only. |
| Isotype | IgG |
| Clonality | Polyclonal |
| Purity | Antigen Affinity Purified |
| Uniprot ID | P01019 |
| Dilution Range | WB:1:500-2000, IHC-P:1:20-1:200 |
| Expiration Date | 12 months from date of receipt. |



Immunofluorescence analysis of HepG2 cells using Angiotensin II antibody