

Product Datasheet

SQSTM1 antibody orb100344

Description	Rabbit Polyclonal to SQSTM1.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Predicted Reactivity	Bovine, Canine, Equine, Porcine, Sheep
Conjugation	Unconjugated
Custom conjugations	Unconjugated, HRP
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived between 61-140 amino acids of human SQSTM1
Target	SQSTM1
Form/Appearance	Liquid: 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Concentration	- 100 µg (in 100 µl): 1 mg/ml - 200 µg (in 200 µl): 1 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified by Protein A
MW	48 kDa
Uniprot ID	Q13501 (Human); Q64337 (Mouse); O08623 (Rat)
Entrez	8878
Dilution Range	WB: 1:500-2000, ELISA: 1:500-1000, IHC-P: 1:100-500, IF: 1:100-500
Expiration Date	12 months from date of receipt.