


Product Datasheet

IgG2b Isotype control (PE-Cy5) orb12200

Description	Mouse monoclonal to IgG2b Isotype control. This can be used as a mouse IgG2b Isotype control in flow cytometry and other applications.
Species/Host	Mouse
Conjugation	PE/Cy5
Tested Applications	FACS, ICC, IF
Target	IgG2b
Form/Appearance	Phosphate-buffered solution, pH 7.4, containing 0.09% sodium azide and 0.2% (w/v) BSA
Storage	Store at 4°C. Do not freeze. Light sensitive material.
Note	For research use only.
Isotype	IgG2b
Clonality	Monoclonal
Purity	Affinity purified
Clone ID	4A2B5
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of peripheral blood lymphocytes using IgG2b Isotype Control antibody (PE-Cy5)