

Product Datasheet

IVD antibody (Biotin) orb450722

Description	Rabbit polyclonal antibody against IVD conjugated to Biotin
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Biotin
Tested Applications	ELISA, ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human IVD
Target	IVD (Biotin)
Form/Appearance	Liquid: 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified by Protein A
MW	43 kDa
Uniprot ID	P26440 (Human); Q9JHI5 (Mouse); P12007 (Rat)
Entrez	3712
Dilution Range	WB: 1:50-200, ELISA: 1:100-1000, IHC-P: 1:50-200, ICC: 1:50-200, IF: 1:50-200
Expiration Date	12 months from date of receipt.