

Product Datasheet

PNR antibody (Biotin) orb453088

Description	Rabbit polyclonal antibody against PNR conjugated to Biotin
Species/Host	Rabbit
Conjugation	Biotin
Tested Applications	ELISA, ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from mouse PNR
Target	PNR (Biotin)
Form/Appearance	Liquid: 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified by Protein A
MW	38 kDa
Uniprot ID	O14804 (Human); Q5QD14 (Mouse); Q5QD23 (Rat)
Entrez	9038
Dilution Range	WB: 1:50-200, ELISA: 1:100-1000, IHC-P: 1:50-200, ICC: 1:50-200, IF: 1:50-200
Expiration Date	12 months from date of receipt.