

Product Datasheet

RSK3 (phospho-Thr353/356) antibody (Biotin) (orb502226)

Description	Rabbit polyclonal antibody against RSK3 (phospho-Thr353/356) conjugated to Biotin...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Biotin
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human RSK3 around the phosphorylation site of Thr353/356
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	83kDa
Uniprot ID	Q9WUT3 , Q15349
Dilution Range	WB: 1:50-200, ELISA: 1:100-1000, IHC-P: 1:50-200, IF: 1:50-200
Expiration Date	12 months from date of receipt.