

Product Datasheet

GLP1 peptide orb99760

Description	GLP1 peptide
Reactivity	Bovine, Human
Conjugation	Unconjugated
Tested Applications	HPLC
Target	GLP1
Storage	Aliquot and store at -20°C. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Protein Sequence	Glp-Arg-Pro-Arg-Leu-Ser-His-Lys-Gly-Pro-Met-Pro-Phe
Purity	Greater than 98 % purified
Formula	C ₆₉ H ₁₀₈ N ₂₁ O ₁₆ S ₁
MW	1533.81 Da
Expiration Date	12 months from date of receipt.