

## Product Datasheet

### Gr-1 antibody orb433622

<b>Description</b>	Rat monoclonal antibody to Gr-1
<b>Species/Host</b>	Rat
<b>Reactivity</b>	Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	FACS, IHC-Fr, IP, WB
<b>Target</b>	Gr-1 antibody
<b>Form/Appearance</b>	Liquid: Phosphate buffered saline
<b>Concentration</b>	1.0 mg/ml
<b>Preservatives</b>	0.09% Sodium Azide
<b>Storage</b>	Store at 4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
<b>Note</b>	For research use only.
<b>Isotype</b>	IgG2b
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Purified IgG
<b>Clone ID</b>	RB6-8C5
<b>Uniprot ID</b>	<b>P35461</b>
<b>Expiration Date</b>	12 months from date of receipt.