

Product Datasheet

Melanostatin peptide orb104168

Description	Melanostatin peptide
Reactivity	Frog
Conjugation	Unconjugated
Tested Applications	HPLC
Target	Melanostatin
Storage	Aliquot and store at -20°C. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Protein Sequence	Tyr-Pro-Ser-Lys-Pro-Asp-Asn-Pro-Gly-Glu-Asp-Ala-Pro-Ala-Glu-Asp-Met-Ala-Lys-Tyr-Tyr-Ser-Ala-Leu-Arg-His-Tyr-Ile-Asn-Leu-Ile-Thr-Arg-Gln-Arg-Tyr-NH ₂
Purity	Greater than 98 % purified
Formula	C ₁₈₉ H ₂₈₅ N ₅₃ O ₅₇ S ₁
MW	4243.8 Da
Expiration Date	12 months from date of receipt.