

Product Datasheet

GluN2B/NR2B antibody (orb148746)

Description	Mouse monoclonal antibody against GRIN2B conjugated to ATTO-390...
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	ATTO-390
Tested Applications	AM, ICC, IP, WB, IHC
Immunogen	Fusion protein amino acids 20-271 (extracellular N-terminus) of rat NR2B
Target	NR2B
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only.
Isotype	IgG2b
Clonality	Monoclonal
MW	166kDa
Clone ID	S59
Uniprot ID	Q00960
Entrez	24410
Dilution Range	WB (1:1000), IHC (1:1000), ICC/IF (1:100)
Expiration Date	12 months from date of receipt.