

## Product Datasheet

### GluN2B/NR2B antibody (orb148749)

<b>Description</b>	Mouse monoclonal antibody against GRIN2B conjugated to ATTO-594...
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	ATTO-594
<b>Tested Applications</b>	AM, ICC, IP, WB, IHC
<b>Immunogen</b>	Fusion protein amino acids 20-271 (extracellular N-terminus) of rat NR2B
<b>Target</b>	NR2B
<b>Concentration</b>	1 mg/ml
<b>Storage</b>	Conjugated antibodies should be stored at 4°C
<b>Note</b>	For research use only.
<b>Isotype</b>	IgG2b
<b>Clonality</b>	Monoclonal
<b>MW</b>	166kDa
<b>Clone ID</b>	S59
<b>Uniprot ID</b>	<a href="#">Q00960</a>
<b>Entrez</b>	<a href="#">24410</a>
<b>Dilution Range</b>	WB (1:1000), IHC (1:1000), ICC/IF (1:100)
<b>Expiration Date</b>	12 months from date of receipt.