

Product Datasheet

BRF1 peptide orb217381

Description	The peptide is used to block anti-BRF1 antibody reactivity.
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Target	BRF1
Form/Appearance	Lyophilized powder
Storage	Store at 4°C for short term and at -20°C for long term storage
Note	For research use only.
Purity	>85%
Source	Synthetic
Uniprot ID	Q92994 (Human); Q8CFK2 (Mouse)
Expiration Date	12 months from date of receipt.
Application Notes	This peptide was used as immunogen to produce orb214025 .