

Product Datasheet

COL18A1 peptide orb218660

Description	The peptide is used to block anti-Collagen 18 alpha 1 antibody reactivity.
Reactivity	Canine, Human, Rat
Conjugation	Unconjugated
Target	COL18A1
Form/Appearance	Lyophilized powder
Storage	Store at 4°C for short term and at -20°C for long term storage
Note	For research use only.
Purity	>85%
Source	Synthetic
Uniprot ID	P39060 (Human)
Expiration Date	12 months from date of receipt.
Application Notes	This peptide was used as immunogen to produce orb215304 .