

Product Datasheet

SNRK antibody (PE) orb489960

Description	Rabbit polyclonal antibody against SNRK conjugated to PE
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	PE
Immunogen	KLH conjugated synthetic peptide derived from human SNRK
Target	SNRK (PE)
Form/Appearance	Liquid: 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified by Protein A
MW	84 kDa
Uniprot ID	Q9NRH2 (Human); Q8VDU5 (Mouse); Q63553 (Rat)
Entrez	54861
Expiration Date	12 months from date of receipt.