

Product Datasheet

S100A10 peptide orb217907

Description	S100A10 peptide
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Target	S100A10
Form/Appearance	Lyophilized powder
Storage	Store at 4°C for short term and at -20°C for long term storage
Note	For research use only.
Purity	>85%
Source	Synthetic
Uniprot ID	P60903 (Human); P08207 (Mouse); P05943 (Rat)
Expiration Date	12 months from date of receipt.
Application Notes	This peptide was used as immunogen to produce orb214551 .