

## Product Datasheet

### Collagen I peptide orb217113

<b>Description</b>	Collagen I peptide
<b>Reactivity</b>	Canine, Human, Mouse, Rat, Sheep
<b>Conjugation</b>	Unconjugated
<b>Target</b>	Collagen I
<b>Form/Appearance</b>	Lyophilized powder
<b>Storage</b>	Store at 4°C for short term and at -20°C for long term storage
<b>Note</b>	For research use only.
<b>Purity</b>	>85%
<b>Source</b>	Synthetic
<b>Uniprot ID</b>	<b>P08123</b> (Human); <b>Q01149</b> (Mouse); <b>P02466</b> (Rat)
<b>Expiration Date</b>	12 months from date of receipt.
<b>Application Notes</b>	This peptide was used as immunogen to produce <b>orb213757</b> .