

Product Datasheet

CREB1 peptide orb217131

Description	CREB1 peptide
Reactivity	Feline, Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Target	CREB1
Form/Appearance	Lyophilized powder
Storage	Store at 4°C for short term and at -20°C for long term storage
Note	For research use only.
Purity	>85%
Source	Synthetic
Uniprot ID	P16220 (Human); Q01147 (Mouse); P15337 (Rat)
Expiration Date	12 months from date of receipt.
Application Notes	This peptide was used as immunogen to produce orb213775 .