

Product Datasheet

CTNNB1 peptide orb217157

Description	The peptide is used to block anti-Beta-catenin antibody reactivity.
Reactivity	Canine, Feline, Human, Monkey, Mouse, Rabbit, Rat, Sheep, Zebrafish
Conjugation	Unconjugated
Target	CTNNB1
Form/Appearance	Lyophilized powder
Storage	Store at 4°C for short term and at -20°C for long term storage
Note	For research use only.
Purity	>85%
Source	Synthetic
Uniprot ID	P35222 (Human); Q02248 (Mouse); Q9WU82 (Rat)
Expiration Date	12 months from date of receipt.
Application Notes	This peptide was used as immunogen to produce orb213801 .