

Product Datasheet

COPS3 peptide orb218157

Description	The peptide is used to block anti-COPS3 antibody reactivity.
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Target	COPS3
Form/Appearance	Lyophilized powder
Storage	Store at 4°C for short term and at -20°C for long term storage
Note	For research use only.
Purity	>85%
Source	Synthetic
Uniprot ID	Q9UNS2 (Human); O88543 (Mouse); Q68FW9 (Rat)
Expiration Date	12 months from date of receipt.
Application Notes	This peptide was used as immunogen to produce orb214801 .