

Product Datasheet

Calprotectin antibody orb108131

Description	Mouse monoclonal antibody to Calprotectin
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Form/Appearance	Lyophilized
Preservatives	NaN3
Note	For research use only.
Clonality	Monoclonal
Purity	Purified
Source	Cell Culture
Expiration Date	32 months from date of receipt.
Application Notes	Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.