

Product Datasheet

14-3-3 antibody (HRP) orb109986

Description	Rabbit polyclonal to 14-3-3 (HRP).
Species/Host	Rabbit
Conjugation	HRP
Tested Applications	ELISA, IHC-P, WB
Immunogen	Synthetic peptide (Please contact us for more details)
Target	14-3-3
Form/Appearance	Liquid: 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified by Protein A
MW	30 kDa
Dilution Range	WB: 1:50-200, ELISA: 1:100-1000, IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.
Application Notes	orb109986 is the HRP conjugated form of orb10003