


Product Datasheet

Human, Bovine, Porcine, Rat [Ala11, 22, 28] VIP peptide orb233712

| | |
|------------------------|---|
| Description | Human, Bovine, Porcine, Rat [Ala11, 22, 28] VIP peptide |
| Reactivity | Human, Rat |
| Conjugation | Unconjugated |
| Sequence | His-Ser-Asp-Ala-Val-Phe-Thr-Asp-Asn-Tyr-Ala-Arg-Leu-Arg-Lys-Gln-Met-Ala-Val-Lys-Lys-Ala-Leu-Asn-Ser-Ile-Leu-Ala-NH ₂ |
| Target | [Ala11, 22, 28] VIP |
| Form/Appearance | Lyophilized powder |
| Storage | Store at 4°C for short term and at -20°C for long term storage |
| Note | For research use only. |
| Purity | > 95% |
| Formula | C ₁₃₉ H ₂₃₁ N ₄₃ O ₃₉ S |
| MW | 3160.72 |
| Source | Synthetic |
| Expiration Date | 12 months from date of receipt. |



Human, Bovine, Porcine,
Rat
[Ala11, 22, 28] VIP peptide