


Product Datasheet

beta-Amyloid-Protein [Arg15, Asp16, 25, Pro18, 21, 23, Val22, Ile24] peptide orb233889

Description	beta-Amyloid-Protein [Arg15, Asp16, 25, Pro18, 21, 23, Val22, Ile24] peptide
Conjugation	Unconjugated
Sequence	Arg-Asp-Leu-Pro-Phe-Phe-Pro-Val-Pro-Ile-Asp
Target	beta-Amyloid-Protein [Arg15, Asp16, 25, Pro18, 21, 23, Val22, Ile24]
Form/Appearance	Lyophilized powder
Storage	Store at 4°C for short term and at -20°C for long term storage
Note	For research use only.
Purity	> 95%
Formula	C ₆₄ H ₉₄ N ₁₄ O ₁₆
MW	1315.55
Source	Synthetic
Expiration Date	12 months from date of receipt.



beta-Amyloid-Protein
[Arg15,
Asp16, 25, Pro18, 21, 23,
Val22, Ile24] peptide