

Product Datasheet

Calbindin antibody orb389562

Description	Mouse monoclonal antibody to Calbindin
Species/Host	Mouse
Reactivity	Bovine, Equine, Gallus, Human, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Immunogen	Full-length recombinant human protein
Target	Calbindin
Form/Appearance	Liquid
Concentration	1 mg/ml
Storage	Aliquoted and stored frozen at -20°C to -70°C in a manual defrost freezer for six months without detectable loss of activity. The antibody can be stored at 2°-8°C for 1 month without detectable loss of activity. Avoid repeated freeze-thaw cycles.
Note	For research use only.
Isotype	IgG1
Clonality	Monoclonal
Purity	Affinity Purified
Dilution Range	IF/ICC: 1:50000, IHC-P: 1:5000, WB: 1:5000
Expiration Date	12 months from date of receipt.