

Product Datasheet

CLN3 antibody orb12725

Description	Rabbit polyclonal antibody to CLN3
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant protein of human CLN3
Target	CLN3
Form/Appearance	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	48 kDa
Uniprot ID	Q13286
Entrez	1201
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:100
Expiration Date	12 months from date of receipt.

