

Product Datasheet

ACE2 peptide orb339324

Description	The peptide is used to ACE2 antibody reactivity.
Conjugation	Unconjugated
Target	ACE2
Form/Appearance	Lyophilized powder
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Source	Synthetic
Expiration Date	12 months from date of receipt.
Application Notes	This peptide was used as immunogen to produce orb338930 (ACE2 antibody) The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.