

Product Datasheet

14-3-3 antibody (PE) orb124427

Description	Rabbit polyclonal to 14-3-3 (PE).
Species/Host	Rabbit
Conjugation	PE
Tested Applications	IF
Immunogen	KLH conjugated synthetic peptide derived from Malus domestica 14-3-3 family protein C-terminus
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	30kDa
Dilution Range	IF: 1:50-200
Expiration Date	12 months from date of receipt.
Application Notes	orb124427 is the PE conjugated form of orb10003