

Product Datasheet

CD133 peptide orb374954

Description	CD133 peptide
Conjugation	Unconjugated
Tested Applications	HPLC
Target	CD133
Form/Appearance	Lyophilized
Storage	Upon reconstitution, aliquot and store at -20°C or -80°C. for long term storage, Avoid repeated freeze/thaw cycles.
Note	For research use only.
Source	Synthesis
Expiration Date	12 months from date of receipt.
Application Notes	This peptide was used as immunogen to produce orb99113 (CD133 antibody)