



## Product Datasheet

### AMBRA1 antibody orb96241

<b>Description</b>	Rabbit Polyclonal to AMBRA1.
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Bovine, Gallus, Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC-P, WB
<b>Target</b>	AMBRA1
<b>Form/Appearance</b>	PBS,0.05% BSA,0.05% sodium azide
<b>Concentration</b>	0.5 mg/ml
<b>Storage</b>	Store at 4°C.stable for 6 months.For long term storage,aliquot and store at -20°C.
<b>Note</b>	For research use only.
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Purity</b>	Protein G Column
<b>Dilution Range</b>	IHC-P:10 µg/ml,WB:0.5-2 µg/ml
<b>Expiration Date</b>	12 months from date of receipt.



Immunohistochemical staining of paraffin embedded human small intestine tissue using AMBRA1 antibody (primary antibody at 1:200)



Western blot analysis of human brain lysate using AMBRA1 antibody