

Product Datasheet

CNOT7 peptide orb382252

Description	The peptide is used to block anti-CNOT7 antibody reactivity
Reactivity	Bovine, Gallus, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	Blocking
Target	CNOT7
Form/Appearance	Lyophilized powder
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Purity	>85%
Source	Synthetic
Uniprot ID	Q9UIV1 (Human); Q60809 (Mouse)
Entrez	29883
Expiration Date	12 months from date of receipt.
Application Notes	This peptide was used as immunogen to produce orb382535 (CNOT7 antibody) Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.