

Product Datasheet

TRIM16L antibody (Biotin) orb455530

Description	Rabbit polyclonal antibody against TRIM16L conjugated to Biotin
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA, ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human TRIM16L
Target	TRIM16L (Biotin)
Form/Appearance	Liquid: 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified by Protein A
MW	40 kDa
Uniprot ID	Q309B1 (Human)
Entrez	147166
Dilution Range	WB: 1:50-200, ELISA: 1:100-1000, IHC-P: 1:50-200, ICC: 1:50-200, IF: 1:50-200
Expiration Date	12 months from date of receipt.