

Product Datasheet

SLC25A32 antibody orb37248

Description	Rabbit polyclonal antibody to SLC25A32
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Custom conjugations	Unconjugated, Biotin, FITC, HRP, CF647, CF488A, CF594, CF405M
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Mitochondrial folate transporter/carrier protein (44-88AA)
Target	SLC25A32
Form/Appearance	Liquid: Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	>95%, Protein G purified
Uniprot ID	Q9H2D1
Dilution Range	WB:1:500-5000
Expiration Date	12 months from date of receipt.

