



Product Datasheet

PHPT1 antibody orb191461

Description	Rabbit polyclonal antibody to PHPT1
Species/Host	Rabbit
Reactivity	Human
Predicted Reactivity	Mouse,Rat,Human
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Immunogen	Recombinant fusion protein of human PHPT1 (NP_054891.2).
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	13kDa
Uniprot ID	Q9NRX4
Entrez	29085
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF/ICC: 1:20-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of various cell lines using PHPT1 antibody.



Confocal immunofluorescence analysis of HeLa cells using PHPT1 antibody.