


Product Datasheet

11A1 (Cytochrome P450) antibody orb229690

Description	Rabbit polyclonal antibody to 11A1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthesized peptide
Target	11A1. Cytochrome P450 11A1 Antibody detects endogenous levels of total Cytochrome P450 11A1
Form/Appearance	PBS, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol
Concentration	1 mg/ml
Storage	Store at 4°C for short term and at -20°C for long term storage
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Immunogen affinity purified
MW	60 kDa
Uniprot ID	P05108
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa cells and A549 cells using 11A1 (Cytochrome P450) antibody