

## Product Datasheet

# Biogenesis of lysosome-related organelles complex 1 subunit 1 antibody (FITC) orb240225

<b>Description</b>	Rabbit polyclonal antibody against Biogenesis of lysosome-related organelles complex 1 subunit 1 conjugated to FITC
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	FITC
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Recombinant Human Biogenesis of lysosome-related organelles complex 1 subunit 1 protein (1-153AA)
<b>Target</b>	BLOC1S1
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Note</b>	For research use only.
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<b>P78537</b>
<b>Dilution Range</b>	ELISA
<b>Expiration Date</b>	12 months from date of receipt.