

Product Datasheet

ZIK1 antibody (FITC) orb222652

| | |
|-----------------------------|--|
| Description | Rabbit polyclonal antibody against ZIK1 conjugated to FITC. |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Predicted Reactivity | Bovine, Canine, Equine, Feline, Porcine |
| Conjugation | FITC |
| Tested Applications | ICC, IF |
| Immunogen | KLH conjugated synthetic peptide derived from human ZIK1 |
| Target | ZIK1 (FITC) |
| Form/Appearance | Liquid: 0.01M TBS (pH7.4), 0.5% BSA, 0.015% ProClin™ 300 and 25% glycerol |
| Concentration | - 100 µg (in 200 µl): 0.5 mg/ml |
| Storage | Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles. |
| Note | For research use only. |
| Isotype | IgG |
| Clonality | Polyclonal |
| Purity | Affinity purified by Protein A |
| MW | 55 kDa |
| Uniprot ID | Q3SY52 (Human); Q80YP6 (Mouse) |
| Entrez | 284307 |
| Dilution Range | ICC: 1:50-200, IF: 1:50-200 |
| Expiration Date | 12 months from date of receipt. |