

## Product Datasheet

### BDNF antibody (FITC) orb222658

<b>Description</b>	Mouse monoclonal antibody against BDNF conjugated to FITC
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human, Mouse, Rat
<b>Predicted Reactivity</b>	Bovine, Canine, Equine, Porcine
<b>Conjugation</b>	FITC
<b>Tested Applications</b>	ICC, IF
<b>Immunogen</b>	Recombined human BDNF full length protein
<b>Target</b>	BDNF (FITC)
<b>Form/Appearance</b>	Liquid: 0.01M TBS (pH7.4), 0.5% BSA, 0.015% ProClin™ 300 and 25% glycerol
<b>Concentration</b>	- 100 µg (in 200 µl): 0.5 mg/ml
<b>Storage</b>	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
<b>Note</b>	For research use only.
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Affinity purified by Protein G
<b>MW</b>	13/27 kDa
<b>Clone ID</b>	3B2
<b>Uniprot ID</b>	<a href="#">P23560</a> (Human); <a href="#">P21237</a> (Mouse); <a href="#">P23363</a> (Rat)
<b>Entrez</b>	<a href="#">627</a>
<b>Dilution Range</b>	ICC: 1:50-200, IF: 1:50-200
<b>Expiration Date</b>	12 months from date of receipt.