

## Product Datasheet

### IVD antibody (HRP) orb517219

<b>Description</b>	Rabbit polyclonal antibody against IVD conjugated to HRP
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	HRP
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Recombinant Human Isovaleryl-CoA dehydrogenase, mitochondrial protein (139-241AA)
<b>Target</b>	Recognizes endogenous levels of IVD
<b>Form/Appearance</b>	Liquid: Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Concentration</b>	10 mM - 11 mM
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Note</b>	For research use only.
<b>Reconstitution</b>	Reconstitute with 200 ul of sterile H <sub>2</sub> O.
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Purity</b>	>95%, Protein G purified
<b>MW</b>	878.67 g/mol (free acid)
<b>Uniprot ID</b>	<a href="#">P26440</a>
<b>NCBI</b>	<a href="#">NP_000454.1</a> ; <a href="#">NM_000463.2</a>
<b>Entrez</b>	<a href="#">54658</a>
<b>Dilution Range</b>	WB:1:500-1:5000, IHC-P:1:10-1:100, IF:1:1-1:10
<b>Expiration Date</b>	12 months from date of receipt.