

Product Datasheet

IVD antibody (FITC) orb517220

Description	Rabbit polyclonal antibody against IVD conjugated to FITC
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Tested Applications	ELISA
Immunogen	Recombinant Human Isovaleryl-CoA dehydrogenase, mitochondrial protein (139-241AA)
Target	Recognizes endogenous levels of IVD
Form/Appearance	Liquid: Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Concentration	10 mM - 11 mM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only.
Reconstitution	Reconstitute with 200 ul of sterile H2O.
Isotype	IgG
Clonality	Polyclonal
Purity	>95%, Protein G purified
MW	878.67 g/mol (free acid)
Uniprot ID	P26440
NCBI	NP_000454.1; NM_000463.2
Entrez	54658
Dilution Range	WB:1:500-1:5000, IHC-P:1:10-1:100, IF:1:1-1:10
Expiration Date	12 months from date of receipt.