

Product Datasheet

Neuropeptide Y antibody (FITC) orb559281

Description	Rabbit polyclonal antivosy against Neuropeptide Y conjugated to FITC
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Tested Applications	IF
Immunogen	KLH conjugated synthetic peptide derived from human C-flanking Peptide Of NPY
Target	Neuropeptide Y
Form/Appearance	Liquid: 0.01M TBS (pH7.4), 0.5% BSA, 0.015% ProClin™ 300 and 25% glycerol
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified by Protein A
MW	3.3/8 kDa
Uniprot ID	P01303 (Human); P57774 (Mouse); P07808 (Rat)
Entrez	4852
Dilution Range	IF: 1:50-200
Expiration Date	12 months from date of receipt.