

Product Datasheet

Thymidine Kinase 1 antibody (PE) orb559817

Description	Rabbit polyclonal antivosy against Thymidine Kinase 1 conjugated to PE
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	PE
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human Thymidine Kinase 1
Target	Thymidine Kinase 1
Form/Appearance	Liquid: 0.01M TBS (pH7.4), 0.5% BSA, 0.015% ProClin™ 300 and 25% glycerol
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified by Protein A
MW	25 kDa
Uniprot ID	P04183 (Human)
Entrez	7083
Dilution Range	ICC: 1:50-200, IF: 1:50-200
Expiration Date	12 months from date of receipt.