

Product Datasheet

ZNF608 antibody (PE) orb559858

Description	Rabbit polyclonal antivosy against ZNF608 conjugated to PE
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	PE
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human ZNF608
Target	ZNF608
Form/Appearance	Liquid: 0.01M TBS (pH7.4), 0.5% BSA, 0.015% ProClin™ 300 and 25% glycerol
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified by Protein A
MW	162 kDa
Uniprot ID	Q9ULD9 (Human); Q56A10 (Mouse)
Entrez	57507
Dilution Range	ICC: 1:50-200, IF: 1:50-200
Expiration Date	12 months from date of receipt.