

Product Datasheet

CD41 antibody (PE) orb559898

| | |
|----------------------------|--|
| Description | Mouse monoclonal antibody against CD41 conjugated to PE |
| Species/Host | Mouse |
| Reactivity | Human |
| Conjugation | PE |
| Tested Applications | ICC, IF |
| Immunogen | KLH conjugated synthetic peptide derived from human CD41 |
| Target | CD41 |
| Form/Appearance | Liquid: 0.01M TBS (pH7.4), 0.5% BSA, 0.015% ProClin™ 300 and 25% glycerol |
| Concentration | - 100 µg (in 200 µl): 0.5 mg/ml |
| Storage | Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles. |
| Note | For research use only. |
| Isotype | IgG |
| Clonality | Monoclonal |
| Purity | Affinity purified by Protein G |
| MW | 125 kDa |
| Clone ID | 9B10 |
| Uniprot ID | P08514 (Human); Q9QUM0 (Mouse) |
| Entrez | 3674 |
| Dilution Range | ICC: 1:50-200, IF: 1:50-200 |
| Expiration Date | 12 months from date of receipt. |