

Product Datasheet

RECK antibody (PE) orb559921

Description	Rabbit polyclonal antivosy against RECK conjugated to PE
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	PE
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human RECK
Target	RECK
Form/Appearance	Liquid: 0.01M TBS (pH7.4), 0.5% BSA, 0.015% ProClin™ 300 and 25% glycerol
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified by Protein A
MW	107 kDa
Uniprot ID	Q95980 (Human); Q5W0K6 (Human)
Entrez	313488
Dilution Range	ICC: 1:50-200, IF: 1:50-200
Expiration Date	12 months from date of receipt.