

Product Datasheet

STAT3 antibody (PE) orb559923

Description	Mouse monoclonal antibody against STAT3 conjugated to PE
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Predicted Reactivity	Mouse
Conjugation	PE
Tested Applications	ICC, IF
Immunogen	Recombinant human STAT3 Protein
Target	STAT3
Form/Appearance	Liquid: 0.01M TBS (pH7.4), 0.5% BSA, 0.015% ProClin™ 300 and 25% glycerol
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Monoclonal
Purity	Affinity purified by Protein G
MW	88 kDa
Clone ID	5H4
Uniprot ID	P40763 (Human); P42227 (Mouse); P52631 (Rat)
Dilution Range	ICC: 1:50-200, IF: 1:50-200
Expiration Date	12 months from date of receipt.