

Product Datasheet

FoxO3a (phospho-Ser294) antibody (PE) orb560058

Description	Rabbit polyclonal antivosy against FoxO3a (phospho-Ser294) conjugated to PE
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	PE
Tested Applications	IF
Immunogen	KLH conjugated synthetic peptide derived from the region around the phosphorylation site of S human FoxO3a (phospho-Ser294)
Target	FoxO3a (phospho-Ser294)
Form/Appearance	Liquid: 0.01M TBS (pH7.4), 0.5% BSA, 0.015% ProClin™ 300 and 25% glycerol
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified by Protein A
MW	71 kDa
Uniprot ID	O43524 (Human); Q9WVH4 (Mouse)
Entrez	2309
Dilution Range	IF: 1:50-200
Expiration Date	12 months from date of receipt.