

Product Datasheet

BRLF1 antibody (Biotin) **orb560290**

Description	Rabbit polyclonal antivosy against BRLF1 conjugated to Biotin
Species/Host	Rabbit
Conjugation	Biotin
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human BRLF1
Target	BRLF1
Form/Appearance	Liquid: 0.01M TBS (pH7.4), 0.5% BSA, 0.015% ProClin™ 300 and 25% glycerol
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified by Protein A
MW	66.5 kDa
Dilution Range	WB: 1:50-200, ELISA: 1:100-1000, IHC-P: 1:50-200, IF: 1:50-200
Expiration Date	12 months from date of receipt.