

## Product Datasheet

### CD41 antibody (Biotin) orb560492

<b>Description</b>	Mouse monoclonal antibody against CD41 conjugated to Biotin
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Biotin
<b>Tested Applications</b>	ICC, IF, IHC-P
<b>Immunogen</b>	KLH conjugated synthetic peptide derived from human CD41
<b>Target</b>	CD41
<b>Form/Appearance</b>	Liquid: 0.01M TBS (pH7.4), 0.5% BSA, 0.015% ProClin™ 300 and 25% glycerol
<b>Concentration</b>	- 100 µg (in 200 µl): 0.5 mg/ml
<b>Storage</b>	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
<b>Note</b>	For research use only.
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Affinity purified by Protein G
<b>MW</b>	125 kDa
<b>Clone ID</b>	9B10
<b>Uniprot ID</b>	<b>P08514</b> (Human); <b>Q9QUM0</b> (Mouse)
<b>Entrez</b>	<b>3674</b>
<b>Dilution Range</b>	IHC-P: 1:50-200, ICC: 1:50-200, IF: 1:50-200
<b>Expiration Date</b>	12 months from date of receipt.