

Product Datasheet

Beta Catenin antibody (Biotin) orb560526

Description	Mouse monoclonal antibody against Beta Catenin conjugated to Biotin
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Biotin
Tested Applications	ICC, IF, IHC-P, WB
Immunogen	Recombinant human β -Catenin
Target	Beta Catenin
Form/Appearance	Liquid: 0.01M TBS (pH7.4), 0.5% BSA, 0.015% ProClin™ 300 and 25% glycerol
Concentration	- 100 μ g (in 200 μ l): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Monoclonal
Purity	Affinity purified by Protein G
MW	92 kDa
Clone ID	5D7
Uniprot ID	P35222 (Human); Q02248 (Mouse); Q9WU82 (Rat)
Entrez	1499
Dilution Range	WB: 1:50-200, IHC-P: 1:50-200, ICC: 1:50-200, IF: 1:50-200
Expiration Date	12 months from date of receipt.