

Product Datasheet

Angiotensin antibody (Biotin) orb560600

Description	Rabbit polyclonal antiserum against Angiotensin conjugated to Biotin
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Biotin
Tested Applications	ELISA, ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Angiotensin
Target	Angiotensin
Form/Appearance	Liquid: 0.01M TBS (pH7.4), 0.5% BSA, 0.015% ProClin™ 300 and 25% glycerol
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified by Protein A
MW	118 kDa
Uniprot ID	Q4VCS5 (Human); Q8VHG2 (Mouse)
Entrez	154796
Dilution Range	WB: 1:50-200, ELISA: 1:100-1000, IHC-P: 1:50-200, ICC: 1:50-200, IF: 1:50-200
Expiration Date	12 months from date of receipt.